


Product Datasheet

GRIK2 antibody (orb676223)

Description	GRIK2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic Peptide
Target	GRIK2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13002
Dilution Range	IHC-P: 1: 100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of Human
Co...