



Product Datasheet

KCNH5 antibody (orb676232)


Description	KCNH5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic Peptide
Target	KCNH5
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NCM2
Dilution Range	WB: 1: 1000-2000, IHC-P: 1: 100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse Brain



Western blot analysis of 1) Rat Brain Tissue



Immunohistochemical staining of Rat Brain