



Product Datasheet

MAPK14 antibody (orb676236)


Description	MAPK14 antibody
Species/Host	Rabbit
Reactivity	Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the C-terminal region of Human p38.
Target	MAPK14
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16539
Dilution Range	WB: 1:500-1:2000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse-brain lys...



Western blot analysis of HELA KB PC-3 us...



Western blot analysis of HELA HEPG2 PC-3...