



Product Datasheet

F13B antibody (orb676258)

Description	F13B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from the Internal region of Human Factor XIII B.
Target	F13B
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05160
Dilution Range	WB: 1:500-1:2000, IHC-P: 1: 100-1: 300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of various lysates us...



Immunohistochemical staining of human-li...