



Product Datasheet

ACER2 antibody (orb676276)


Description	ACER2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from the Internal region of Human ASAH3L. at AA range: 50-130
Target	ACER2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5QJU3
Dilution Range	WB: 1:500-1:2000, IHC-P: 1: 100-1: 300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human-li...



Western blot analysis of COS7 using ACER...



Immunohistochemical staining of human-br...