



Product Datasheet

MET antibody (orb676285)


Description	MET antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Protein of Met of MET
Target	MET
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08581
Dilution Range	IHC-P: 1: 100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Liver tissue



Immunohistochemical staining of Human Lung tissue



Immunohistochemical staining of Human Stomach tissue