


Product Datasheet

PMEPA1 antibody (orb676358)

Description	PMEPA1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide from Human protein at AA range: 1268-1341
Target	PMEPA1
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q969W9
Dilution Range	WB: 1: 500-1: 2000, ELISA: 1: 5000-1: 20000
Expiration Date	12 months from date of receipt.



Western blot analysis of 823 22RV1 HELA ...