



## Product Datasheet

### FCGR3A antibody (orb676390)


|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| <b>Description</b>         | FCGR3A antibody                                                              |
| <b>Species/Host</b>        | Rabbit                                                                       |
| <b>Reactivity</b>          | Human                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                 |
| <b>Tested Applications</b> | ELISA, IHC                                                                   |
| <b>Immunogen</b>           | Synthetic peptide from Human protein at AA range: 100-150                    |
| <b>Target</b>              | FCGR3A/FCGR3B                                                                |
| <b>Preservatives</b>       | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| <b>Form/Appearance</b>     | Liquid                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                |
| <b>Note</b>                | For research use only                                                        |
| <b>Isotype</b>             | IgG                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                   |
| <b>Uniprot ID</b>          | <b>O75015</b> , <b>P08637</b>                                                |
| <b>Dilution Range</b>      | IHC-P: 1: 50-300, ELISA: 1: 10000-20000                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                              |



Immunohistochemical staining of human-br...



Immunohistochemical staining of human-br...



Immunohistochemical staining of human-ki...