


Product Datasheet

WIPI2 antibody (orb676514)

Description	WIPI2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Protein
Target	WIPI2
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y4P8
Dilution Range	IHC-P: 1:50-1:200, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Br...



Immunohistochemical staining of Human Sk...