



Product Datasheet

HSPA5 antibody (orb676526)


Description	HSPA5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide from Human protein at AA range: 1268-1341
Target	HSPA5
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11021
Dilution Range	WB: 1: 500-1: 2000, IHC-P: 1: 50-1: 300, ELISA: 1: 5000-1: 20000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human-
co...



Western blot analysis of 3T3 Cell lysate...



Immunohistochemical staining of human-
br...