

## Product Datasheet

### GUCY2D antibody (orb678679)

|                            |                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------|
| <b>Description</b>         | GUCY2D antibody                                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                          |
| <b>Tested Applications</b> | ELISA, WB                                                                             |
| <b>Immunogen</b>           | Human GUCY2D                                                                          |
| <b>Target</b>              | GUCY2D                                                                                |
| <b>Preservatives</b>       | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Form/Appearance</b>     | Liquid                                                                                |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                         |
| <b>Note</b>                | For research use only                                                                 |
| <b>Isotype</b>             | IgG                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                            |
| <b>Uniprot ID</b>          | <a href="#">Q02846</a>                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                       |