



Product Datasheet

Peroxiredoxin 1 antibody (orb6796)


Description	Peroxiredoxin 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human PRDX1 (1-100/199 aa)
Target	Peroxiredoxin 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q06830
Dilution Range	WB=1:500-2000, Flow-Cyt=2µg/Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human Raji cell; Lane...



Flow cytometry of mouse splenocytes (blue...)



Flow cytometry of HeLa. Primary Antibody ...