

Product Datasheet

PPP1R9B antibody (orb680460)

Description	PPP1R9B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human PPP1R9B
Target	PPP1R9B
Preservatives	0.1M NaHCO ₃ , 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96SB3
Expiration Date	12 months from date of receipt.