

Product Datasheet

RAB31 antibody (orb680679)

Description	RAB31 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human RAB31
Target	RAB31
Preservatives	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13636
Dilution Range	WB:1:500-1:1000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.