

## Product Datasheet

### XYLT1 antibody (orb682302)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | XYLT1 antibody                                                                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                         |
| <b>Tested Applications</b> | ELISA, WB                                                                                                            |
| <b>Immunogen</b>           | Human XYLT1                                                                                                          |
| <b>Target</b>              | XYLT1                                                                                                                |
| <b>Preservatives</b>       | 0.1M NaHCO <sub>3</sub> , 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Form/Appearance</b>     | Liquid                                                                                                               |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                        |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Isotype</b>             | IgG                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                           |
| <b>Uniprot ID</b>          | <b>Q86Y38</b>                                                                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |