

Product Datasheet

PDE9A antibody (orb682790)

Description	PDE9A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human PDE9A protein
Target	PDE9A
Preservatives	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	076083
Expiration Date	12 months from date of receipt.