


## Product Datasheet

### MSR1 antibody (orb682975)

|                            |                                                                                                                   |
|----------------------------|-------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | MSR1 antibody                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                      |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                    |
| <b>Immunogen</b>           | Recombinant protein of Human MSR1                                                                                 |
| <b>Target</b>              | MSR1                                                                                                              |
| <b>Preservatives</b>       | Store at -20°C or -80°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3. |
| <b>Form/Appearance</b>     | liquid                                                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                     |
| <b>Note</b>                | For research use only                                                                                             |
| <b>Isotype</b>             | IgG                                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                                        |
| <b>Uniprot ID</b>          | <b>P21757</b>                                                                                                     |
| <b>Dilution Range</b>      | WB: 1: 500-1: 2000, IHC-P: 1: 50-1: 200                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                   |



Western blot analysis of Various cell li...