


Product Datasheet

Cleaved-PARP1 (A214) antibody (orb683081)

Description	Cleaved-PARP1 (A214) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PARP1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 °C.Stable for 14 months from date of receipt
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09874
Dilution Range	WB: 1: 500-1: 2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cells using...