


Product Datasheet

DDIT3 antibody (orb683089)

Description	DDIT3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	A synthesized peptide derived from Human GADD153
Target	DDIT3
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35638
Dilution Range	WB: 1: 500-1: 2000, IHC-P: 1: 50-1: 200, IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...