


Product Datasheet

EWS antibody (orb683216)

Description	EWS antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.213~217(T-Y-G-Q-P)derived from Human EWS.
Target	EWS
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZD1
Dilution Range	WB: 1: 500-1: 1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 and 3T3 cel...