


Product Datasheet

Synpr antibody (orb683220)

Description	Synpr antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.257~261(G-P-T-S-F)derived from Rat synaptophysin 2.
Target	Synpr
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22831
Dilution Range	WB: 1: 500-1: 1000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain and mouse brain.