

Product Datasheet

CUEDC2 antibody (orb683280)

Description	CUEDC2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.262~266(D-V-R-N-P) derived from Human CUEDC2.
Target	CUEDC2
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H467
Dilution Range	WB: 1: 500-1: 1000
Expiration Date	12 months from date of receipt.

