


Product Datasheet

E2F2 antibody (orb683313)

Description	E2F2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.425~429(L-F-D-S-Y) derived from Human E2F-2.
Target	E2F2
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14209
Dilution Range	WB: 1: 500-1: 1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa and K562 c...