


Product Datasheet

GFI1B antibody (orb683315)

Description	GFI1B antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.322~326(H-R-E-S-Q) derived from Human GFI1b.
Target	GFI1B
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5VTD9
Dilution Range	WB: 1: 500-1: 1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji cells using...