


Product Datasheet

MC3R antibody (orb684610)

Description	MC3R antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide corresponding to residues near the N terminal of Human melanocortin 3 receptor
Target	MC3R
Preservatives	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41968
Dilution Range	WB: 1: 500-1: 3000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse large int...