

Product Datasheet

Fbxl4 antibody (orb684618)

Description	Fbxl4 antibody
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around aa.1~5 (M-S-L-L-A) derived from Fruit fly F box and leucine-rich-repeat gene 4.
Target	Fbxl4
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9VY46
Dilution Range	WB: 1: 500-1: 3000
Expiration Date	12 months from date of receipt.