



Product Datasheet

AAK1 antibody (orb684783)

Description	AAK1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fusion protein corresponding to N terminal 236 amino acids of Human Adaptor-associated protein kinase 1
Target	AAK1
Preservatives	Supplied at 1.1mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.3, 0.05% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2M2I8
Dilution Range	WB: 1: 500-1: 1000, IHC-P: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Western blot analysis of Human fetal bra...