



Product Datasheet

MKNK1 (Ab-385) antibody (orb684844)

Description	MKNK1 (Ab-385) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized non-phosphopeptide derived from Human Mnk1 around the phosphorylation site of threonine 385 (L-P-T(p)-P-Q).
Target	MKNK1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BUB5
Dilution Range	WB: 1: 500-1: 3000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells (Laemmli)



Western blot analysis of 3T3 cells using...