



Product Datasheet

SUMO2 antibody (orb685281)

Description	SUMO2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide corresponding to a region derived from internal residues of Human SMT3 suppressor of mif two 3 homolog 2/3/4 (<i>S. cerevisiae</i>)
Target	SUMO2
Preservatives	Rabbit IgG in pH7.3 PBS, 0.05% NaN ₃ , 50% Glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61956
Dilution Range	WB: 1: 500-1: 3000, IHC-P: 1: 50-1: 100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human COX2



Western blot analysis of Human small intestine tissue