


Product Datasheet

Atp1a1 (Ab-23) antibody (orb685383)


Description	Atp1a1 (Ab-23) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Synthesized non-phosphopeptide derived from Human ATP1 α 1/Na+K+ ATPase1 around the phosphorylation site of serine 23 (K-K-S(p)-K-K).
Target	Atp1a1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06685
Dilution Range	WB: 1: 500-1: 3000, IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 cells (Lan...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of HeLa cells usin...