



Product Datasheet

DPYSL2 (Ab-509) antibody (orb685396)


Description	DPYSL2 (Ab-509) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized non-phosphopeptide derived from Human CRMP-2 around the phosphorylation site of threonine 509 (S-V-T(p)-P-K).
Target	DPYSL2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16555
Dilution Range	WB: 1: 500-1: 3000, IHC-P: 1: 50-1: 100
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 cells (Lan...)



Immunohistochemical staining of human br...



Western blot analysis of HT-29 cells usi...