



Product Datasheet

IRS1 (Ab-312) antibody (orb685575)

Description	IRS1 (Ab-312) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Synthesized non-phosphopeptide derived from Human IRS-1 around the phosphorylation site of serine 312 (A-T-S(p)-P-A).
Target	IRS1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35568
Dilution Range	WB: 1: 500-1: 3000, IHC-P: 1: 50-1: 100, IF: 1: 100-1: 500
Expiration Date	12 months from date of receipt.



Western blot analysis of JK cells (Lane ...)



Immunohistochemical staining of human br...



Immunofluorescence analysis of HeLa cell...