



Product Datasheet

TNNI3 antibody (orb685644)


Description	TNNI3 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP
Immunogen	Recombinant Human Troponin I, cardiac muscle protein (2-210AA)
Target	TNNI3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19429
Dilution Range	IHC-P: 1: 20-1: 200, IF: 1: 200-1: 500, IP: 1: 200-1: 2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human he...



Immunoprecipitation analysis of Rat hear...