



Product Datasheet

PTBP3 antibody (orb685745)


Description	PTBP3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Peptide sequence from Human Polypyrimidine tract-binding protein 3 protein (37-55AA)
Target	PTBP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95758
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 20-1: 200, IP: 1: 200-1: 2000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...



Immunohistochemical staining of human ki...



Immunoprecipitation analysis of A549 who...