

Product Datasheet

RTN4RL2 antibody (FITC) (orb685831)

Description	RTN4RL2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Reticulon-4 receptor-like 2 protein (355-372AA)
Target	RTN4RL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86UN3
Expiration Date	12 months from date of receipt.