



Product Datasheet

LRTM2 antibody (orb685849)


Description	LRTM2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Leucine-rich repeat and transmembrane domain-containing protein 2 protein (36-310AA)
Target	LRTM2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N967
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 20-1: 200, IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Western blot analysis of HepG2 whole cel...



Western blot analysis of SH-SY5Y whole c...