

Product Datasheet

VDAC1 antibody (Biotin) (orb685864)

Description	VDAC1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Arabidopsis thaliana Mitochondrial outer membrane protein porin 1 protein (2-276AA)
Target	VDAC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9SRH5
Expiration Date	12 months from date of receipt.