



Product Datasheet

C5AR2 antibody (orb685869)

Description	C5AR2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Peptide sequence from Human C5a anaphylatoxin chemotactic receptor 2 protein (2-19AA)
Target	C5AR2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P296
Dilution Range	IHC-P: 1: 20-1: 200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human placenta



Immunohistochemical staining of human brain