

Product Datasheet

HMGXB4 antibody (Biotin) (orb685876)

Description	HMGXB4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Peptide sequence from Human HMG domain-containing protein 4 protein (489-509AA)
Target	HMGXB4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UGU5
Expiration Date	12 months from date of receipt.