



Product Datasheet

ZRANB1 antibody (orb685877)


Description	ZRANB1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Ubiquitin thioesterase ZRANB1 protein (111-397AA)
Target	ZRANB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UGI0
Dilution Range	WB: 1: 1000-1: 5000, IHC-P: 1: 200-1: 500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human adrenal gland tissue.



Western blot analysis of HeLa whole cell lysates.



Immunohistochemical staining of human prostate tissue.