

Product Datasheet

ART1 antibody (HRP) (orb685898)

Description	ART1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human GPI-linked NAD(P)(+)-arginine ADP-ribosyltransferase 1 protein (46-142AA)
Target	ART1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P52961
Expiration Date	12 months from date of receipt.