

Product Datasheet

ATP2C2 antibody (FITC) (orb686011)

Description	ATP2C2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Calcium-transporting ATPase type 2C member 2 protein (621-702AA)
Target	ATP2C2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075185
Expiration Date	12 months from date of receipt.