

Product Datasheet

Rb (phospho-Ser807/Ser811) antibody (orb6861)

Description	Rb (phospho-Ser807/Ser811) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Rb around the phosphorylation site of Ser807/Ser811 (YI(p-S)PLK(p-S)PY)
Target	Phospho-Rb (Ser807 + Ser811)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102 kDa
Uniprot ID	P06400
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.