

Product Datasheet

At1g32060 antibody (HRP) (orb686193)

Description	At1g32060 antibody (HRP)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Peptide sequence from Arabidopsis thaliana Phosphoribulokinase, chloroplastic protein (223-235AA)
Target	At1g32060
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25697
Expiration Date	12 months from date of receipt.