

Product Datasheet

alaA antibody (FITC) (orb686202)

Description	alaA antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli Glutamate-pyruvate aminotransferase AlaA protein (1-405AA)
Target	alaA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A959
Expiration Date	12 months from date of receipt.