

Product Datasheet

COL9A2 antibody (FITC) (orb686230)

Description	COL9A2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Peptide sequence from Human Collagen alpha-2(IX) chain protein (594-611AA)
Target	COL9A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14055
Expiration Date	12 months from date of receipt.