

Product Datasheet

CARMIL2 antibody (FITC) (orb686254)

Description	CARMIL2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Capping protein, Arp2/3 and myosin-I linker protein 2 protein (1136-1290AA)
Target	CARMIL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6F5E8
Expiration Date	12 months from date of receipt.