

Product Datasheet

SHF antibody (FITC) (orb686274)

Description	SHF antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human SH2 domain-containing adapter protein F protein (1-159AA)
Target	SHF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7M4L6
Expiration Date	12 months from date of receipt.