



Product Datasheet

PARP8 antibody (orb686340)

Description	PARP8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Poly [ADP-ribose] polymerase 8 protein (628-854AA)
Target	PARP8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N3A8
Dilution Range	IHC-P: 1: 200-1: 500, IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon tissue.



Immunofluorescence analysis of HeLa cells.