


Product Datasheet

FCRL5 antibody (orb686364)

Description	FCRL5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Fc receptor-like protein 5 protein (419-553AA)
Target	FCRL5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96RD9
Dilution Range	IF: 1: 50-1: 200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...